## **SPECIFICATIONS:**

Capacity	250-280 samples/hour
Measuring Range	0-15% fat
Accuracy	0.06% for 0-5% fat, 0.10% for 5-8% fat, 0.20% for 8-15% fat
Repeatability	0.03% for 0-5% fat, 0.04% for 5-8% fat, 0.07% for 8-15% fat
Sample Volume	0.5ml/test
Diluent Volume	6.5ml/test
Display type and size	7seg LED, 1/2 " size
Data Interface	BCD OR RS 232. Bluetooth (optional)
Operating temperature	5°C to 45°C
Operating voltage	230V AC +/- 15%, 50Hz or 12V DC
Power Consumption	300w max.
Front keys	Zero, milk in/mix out, repeat, fat
Enclosure	Powder coated aluminium base, MS powder coated top
Dimension	255mm(H) x 460mm(L) X 345MM (W)
Weight	15 Kg